


Introduction

In the development and quality control of biopharmaceuticals (proteins, monoclonal antibodies, antibody drug conjugate, etc.), high performance liquid chromatography (HPLC) is an important tool for analysis and characterisation of their structural heterogeneity. We have developed a novel C4 bonded reversed phase (U)HPLC column, YMC-Triart Bio C4, which is based on organic/inorganic hybrid silica particles with a pore size of 300 Å, designed for biopharmaceuticals separation. Optimised pore size with narrow pore distribution and advanced surface modification that suppresses interaction between an analyte and residual silanol group improve resolution, peak shape, sensitivity and reproducibility of analyses of biomolecules such as intact and subunits of monoclonal antibodies. In this poster, we will show some examples of effective method development for biopharmaceuticals including intact monoclonal antibodies and their fragments with this new hybrid C4 column.


Specifications

Bonded Phase	YMC-Triart Bio C4
USP class	Si-C ₄ H ₉
Base material	inorganic/organic hybrid silica
Particle size (µm)	5, 3, 1.9
Pore size (nm / Å)	30 / 300
Bonding	trifunctional
End-capping	Yes
Usable pH range	1-10
Temperature limit	90°C for pH 1-7 50°C for pH 7-10

Designed for separation of large proteins




Effect of surface modification on peak shape of proteins and peptides



- The chromatograms of protein analysis are compared between YMC-Triart Bio C4 column and its prototype which is different in end-capping process. Although the results obtained with 0.1% TFA (upper figures) are comparable, the results obtained with 0.1% formic acid (lower figures) are significantly different.
- The optimised end-capping process of YMC-Triart Bio C4 suppresses interaction between an analyte and residual silanol group and improves peak shape of proteins with 0.1% formic acid mobile phase.


High chemical stability

Acidic condition (containing 0.1% TFA at 70°C)




- No change in retention time and peak shape is observed when using YMC-Triart Bio C4 under chemically harsh conditions such as the strongly acidic eluents and high temperature.

Comparison of protein separation with 0.1% formic acid mobile phase




Method development for simultaneous analysis of Amyloid β (Aβ) proteins




Analysis of intact monoclonal antibody (mAb)

Effect of increasing temperature




Analysis of five different mAbs




Analysis of mAbs and related substances with LC/MS compatible mobile phase


(A) mAb fragments



(B) Intact mAbs



(C) Antibody-drug conjugate (ADC)



- The analysis results of (A) mAb fragments, (B) intact mAbs, and (C) reduced mAbs and ADCs are compared between 0.1% TFA added condition and 0.1% formic acid added condition. The excellent resolution and peak shape are obtained with YMC-Triart Bio C4 and 0.1% TFA added condition for a variety of mAb and related substances which are different in their molecular size and hydrophobicity.

- Formic acid is a more compatible additive for MS detection than TFA and commonly used in LC/MS analysis of low molecular compounds. However, it usually produces peak broadening and low intensity for proteins as shown in the results with two commercially available C4 columns and (A) mAb fragments. With YMC-Triart Bio C4 column, although slightly broader peaks and shorter retention times are provided using 0.1% formic acid added conditions for larger molecules such as (B) intact mAbs or (C) reduced mAb and ADC, the separation would be very suitable for the structural analysis using LC/MS/(MS).

Conclusions

- The combination of newly developed hybrid particles with a uniform 300 Å pore diameter and advanced surface modification of YMC-Triart Bio C4 column provide excellent peak shape for a variety of proteins and sufficient chemical stability over a wide pH and temperature range. This advantage enables a rapid and efficient method optimisation for complex mixtures of peptides and proteins.
- The superior peak shape and intensity even for larger biopharmaceutical proteins such as intact mAbs, mAb fragments and ADCs, are obtained on YMC-Triart Bio C4 with LC/MS compatible mobile phases containing 0.1% formic acid.