

CHIRAL ART Cellulose-C:
Chiral Separation of the Fungicide Bitertanol




Figure 1: Separation of bitertanol enantiomers on CHIRAL ART Cellulose-C.

Table 1: Chromatographic conditions.

Column	CHIRAL ART Cellulose-C, 3 μ m, 150 x 3 mm ID
Part No.	KCN99S03-1503WT
Eluent	isopropanol / <i>n</i> -hexane / DEA (5/95/0.1)
Flow rate	1 mL/min
Temperature	40°C
Detection	UV at 220 nm
Injection	5 μ L, 1.25 mg/mL solved in isopropanol/ <i>n</i> -hexane (5/95)

Date: 30/01/2017

Author: Dr. Anja Bräutigam