

Baseline separation of chlorophenols as a result of outstanding steric selectivity of YMC-Triart C18 ExRS

Chlorophenols are used as pesticides, bleaching agents and disinfectants. Due to their extreme chemical stability they become concentrated within the food chain. As they are highly toxic their presence in the environment needs to be avoided and monitored

carefully. The various compounds in this group differ only in the amount and/or position of the chlorine substituents. With YMC-Triart C18 ExRS a baseline separation is possible! The high carbon content of this unique U(H)PLC- phase enables steric recognition.




Figure 1: Separation of seven chlorophenols using YMC-Triart C18 ExRS.

Table 1: Method details

| | | |
|--------------------|---|--------------|
| Column | YMC-Triart C18 ExRS 1.9 µm, 8 nm, 75 x 3 mm ID | |
| Part No. | TAR08SP9-L503PT | |
| Eluent | A: Water + 0.1% Formic acid B: Methanol + 0.1% Formic acid | |
| Gradient | Time [min] | Eluent B [%] |
| | 0 | 44 |
| | 8.1 | 50 |
| | 11 | 51.5 |
| | 11.1 | 65 |
| | 20 | 65 |
| Flow rate | 0.7 mL/min | |
| Temperature | 40 °C | |
| Detection | UV at 280 nm | |
| Injection | 1 µL, 0.7 mg/mL each dissolved in methanol | |

Table 2: Chlorophenols analysed

| Retention time | Analyte | Structure |
|----------------|--------------------|---|
| 4.7 min | 4-Chlorophenol |  |
| 7.1 min | 2,6-Dichlorophenol |  |
| 8.9 min | 2,3-Dichlorophenol |  |
| 9.6 min | 2,5-Dichlorophenol |  |
| 10.5 min | 2,4-Dichlorophenol |  |
| 11.1 min | 3,4-Dichlorophenol |  |
| 12.9 min | 3,5-Dichlorophenol |  |