

# Optimisation of oligonucleotide separations using YMC's non porous anion exchanger BioPro IEX QF

Nucleic acid therapeutics such as antisense, siRNA and aptamers are expected to play an important role as the next-generation pharmaceuticals together with antibody drugs.

These drugs demand chromatographic purification





and analysis that can recognize slight structural differences following synthesis.

In this report, useful tips for optimisation of ion exchange chromatography methods for oligonucleotides are provided.

## Samples

1	Single-strand DNA	5'-TCATCACACTGAATACCAAT-3' (DNA 20 mer)
2		5'-GTCATCACACTGAATACCAAT-3' (DNA 21 mer)
3	Single-strand RNA	5'-U(M)C(M)A(M)U(M)C(M)A(M)C(M)A(M)C(M)U(M)G(M)A(M)A(M)U(M)A(M)C(M)C(M)A(M)A(M)U(M)-3' (2'-OMe RNA 20 mer)
4		5'-G(M)U(M)C(M)A(M)U(M)C(M)A(M)C(M)A(M)C(M)U(M)G(M)A(M)A(M)U(M)A(M)C(M)C(M)A(M)A(M)U(M)-3' (2'-OMe RNA 21 mer)
5		5'-UCAUCACACUGAAUACCAAU-3' (RNA 20 mer)
6		5'-GUCAUCACACUGAAUACCAAU-3' (RNA 21 mer)

N(M)=2'-OMe RNA




## Reducing carry over

A) 20 mM Tris-HCl (pH 8.1)  
 B) 20 mM Tris-HCl (pH 8.1) containing 1.0 M NaCl  
 5-70%B (0-15 min), 74%B (15-18 min), 5%B (18-33 min)  
**Initial : 50 mM NaCl**



A) 20 mM Tris-HCl (pH 8.1)  
 B) 20 mM Tris-HCl (pH 8.1) containing 1.0 M NaCl  
 40-70%B (0-15 min), 74%B (15-18 min), 40%B (18-33 min)  
**Initial : 400 mM NaCl**





Column: BioPro IEX QF 5  $\mu$ m, 100 X 4.6 mm ID  
 Part No.: QF00S05-1046WP  
 Flow rate: 1.0 mL/min  
 Temperature: 25  $^{\circ}$ C  
 Detection: UV at 260 nm  
 Injection: 2  $\mu$ L (10 nmol/mL)

Carry over is observed using a gradient with low initial concentration of NaCl. However, good separation with virtually no carry over can be achieved by increasing the initial concentration (e.g. 300-400 mM NaCl).

## Improving peak shape



### 1 Influence of buffer type



Column: BioPro IEX QF (5  $\mu$ m) 100 x 4.6 mm ID  
 Part No.: QF00S05-1046WP  
 Flow rate: 1.0 mL/min  
 Temperature: 25  $^{\circ}$ C  
 Detection: UV at 260 nm  
 Injection: 2  $\mu$ L (10 nmol/mL)  
 Sample: RNA 20 mer

Eluent: A) **X**  
 B) **X** containing 2.0 M NaCl  
 15-100%B (0-30 min)

### 2 Influence of counter ion type






Eluent: A) 10 mM NaOH  
 B) 10 mM NaOH containing 2.0 M **Y**  
 15-100%B (0-30 min) for NaCl  
 5-50%B (0-30 min) for NaClO<sub>4</sub>

Gradient profile is adjusted because eluting strength of NaClO<sub>4</sub> is two to three times more than that of NaCl on ion exchange chromatography.

By changing the buffer from 20 mM Tris-HCl (pH 8.1) to 10 mM NaOH, the tailing factor of the oligonucleotide peak was improved. In addition, changing the counter ion from NaCl to NaClO<sub>4</sub> is effective.

**It is important to optimise buffer and counter ion for excellent peak shape of oligonucleotides.**

## Analysis examples with the optimised conditions



Column: BioPro IEX QF (5 µm) 100 x 4.6 mm ID  
Part No.: QF00S05-1046WP  
Eluent: A) 10 mM NaOH  
B) 10 mM NaOH containing 1.0 M NaClO<sub>4</sub>  
25-55%B (0-15 min), 100%B (15-20 min)  
Flow rate: 1.0 mL/min  
Temperature: 25 °C  
Detection: UV at 260 nm  
Injection: 4 µL (5 nmol/mL each)

Good separation without carry over and peak tailing of oligonucleotides was achieved by optimisation of the buffer/counter ion in the mobile phase and gradient profile, and by using the non porous anion exchange column BioPro IEX QF.