

UHPLC systems compatible with 1 mm ID columns

Fast and high sensitivity UHPLC/MS columns

Ideal choice for


- Low sensitivity compounds
- Low sample amounts
- High sensitivity LC/MS analyses
- Peptides / peptide mapping
- Oligonucleotides

Features

- High performance 1 mm UHPLC columns
- High precision column hardware for reliable results
- Highly reproducible YMC-Triart phases in highly reproducible hardware
- Excellent peak shapes for high sensitivity LC-MS analyses
- Bringing MicroLC and UHPLC advantages together


Specifications	
Particle size	1.9 µm
Pore size	12 nm
Modification	Trifunctional
Base particle	Organic/inorganic hybrid silica
pH range	1–12
Temperature range	pH < 7: 90 °C pH > 7: 50 °C
Pressure limit	100 MPa / 15,000 psi
Recommended flow rate	0.05–0.2 mL/min

YMC-Triart C18



versatile applications
 first choice for
 method development
 pH 1–12/90 °C max.
 100% aqueous eluents

In order to achieve the described benefits, a suitable system for 1 mm UHPLC columns has to be used. Otherwise, the advantages of the columns can be reduced by use of a non-optimised system, with a too high system volume. A selection of suitable UHPLC systems are mentioned below.



Selection of suitable UHPLC systems*

Waters ACQUITY UPLC I-Class (Plus)

For columns up to 4.6 mm ID and 150 mm length

	Sample Manager with Fixed loop (SM-FL)	Sample Manager with flow through needle (SM-FTN)
Total system bandspread (5 σ)	$\leq 7 \mu\text{L}$	$\leq 9 \mu\text{L}$
Dwell volume	95 μL	100 μL
Gradient delay volume	$\leq 75 \mu\text{L}$	
Possible flow rate	0.01–2.0 mL/min	
Maximum pressure	1,200 bar	

Agilent 1290 Infinity (II)

For 1–5 mm ID columns

Has a function to automatically reduce the delay volume

	20–125 μL (depending on configuration)	
Delay volume	Binary pump + sample manager with fixed loop*	20 μL
	Binary pump + Jet Weaver + fixed loop	55 μL
	Binary pump + standard sample manager	90 μL
	Binary pump + Jet Weaver + sample manager	125 μL
Flow rate	0.05–5 mL/min	
Maximum pressure	1,200 bar	

* recommended configuration

Thermo Scientific Vanquish Neo

Uses Viper fittings

1 mm ID columns explicitly mentioned in brochure

Gradient delay volume	$< 2 \mu\text{L}^*$
Flow rate	1 nL/min–0.1 mL/min
Maximum pressure	1,500 bar

* in micro direct injection configuration with 50 μm ID capillaries

Thermo Scientific UltiMate 3000 RSLCnano

System delay volume: $< 350 \text{ nL}$

Uses Viper fittings

	Binary high pressure gradient pump HPG	ternary micro pump TM
Gradient delay volume	$< 25 \text{ nL}$	220 μL
Flow rate	0–0.05 mL/min	0–2.5 mL/min
Maximum pressure	800 bar	620 bar

Selection of suitable microLC systems*

Shimadzu Nexera Mikros

Already used for applications with YMC-Triart C18 1 mm ID columns

For 0.1–1 mm ID columns

Flow rate: 0.1 $\mu\text{L}/\text{min}$ –0.5 mL/min

Uses Universal fittings

SCIEX MicroLC Systems M3 MicroLC, M5 MicroLC, MicroLC 200 Plus

Flow rate: 0.02–0.2 mL/min (different models available with smaller flow rates)

Maximum pressure: 690 bar

Column IDs up to 1 mm possible

*only based on vendor specifications, neither representing a complete selection of systems nor a recommendation for a specific system or vendor.

ACQUITY UPLC and Waters is a trademark of Waters Corp. Agilent is a trademark of Agilent Technologies Inc. Thermo Scientific, Vanquish, UltiMate and Viper are trademarks of Thermo Fisher Scientific Inc. Nexera Mikros is a trademark of Shimadzu Corporation. SCIEX is a trademark of AB SCIEX Pte. Ltd.