



Non-targeted profiling of polar metabolites in human plasma

The fascinating aspect about non-targeted LC analysis is the presence of unknown compounds, for example in biomarker research. A particular challenge for chromatography is then to find the balance between sensitivity and covering a broad spectrum of substances. Hydrophilic liquid chromatography (HILIC) is therefore especially suitable for such analyses

containing polar compounds. It is a perfect match for mass spectrometry (MS) and has very good chromatographic performance. With this combination, non-targeted analysis is well suited to the screening of polar metabolites in human plasma.



This application note shows the non-targeted screening of polar metabolites in human plasma using a bioinert YMC-Accura Triart Diol-HILIC column. The bioinert coated stainless-steel hardware is essential to achieve the highest sensitivity and chromatographic resolution of endogenous isomers.

For this analysis proteins contained in the plasma were precipitated with acetonitrile (3:1 v/v). After decanting and drying the sample, it was re-dissolved with acetonitrile/water (70/30) to the original concentration.

Table 1: Chromatographic conditions [1].

Column:	YMC-Accura Triart Diol-HILIC (1.9 μ m, 12 nm) 150 x 2.1 mm ID
Part No.:	TDH12SP9-15Q1PTC
Eluent:	10 mM ammonium formate in acetonitrile/water (9/1) + 0.1% formic acid 10 mM ammonium formate in acetonitrile/water (1/1) + 0.1% formic acid
Gradient:	1–38%B (0–15 min) [Cleaning step] 99%B (15.1–17 min) [Equilibration] 1%B (17.1–24.9 min)
Flow rate:	0.4 mL/min During equilibration 0.65 mL/min (18–24.8 min)
Temperature:	30 °C
Detection:	ESI positive, Orbitrap Exploris 120, R=60k@mz200 full scan, R=15k DDAtop4
Injection:	1 μ L
Sample:	Human plasma (protein precipitated) in 70% acetonitrile
System:	Thermo Vanquish Horizon



Due to the complex matrix a pre-conditioning step of the stationary phase can be necessary. By using the bioinert YMC-Accura Triart Diol-HILIC column conditioning is already achieved after 4 runs (see Figure 1). Furthermore, an excellent retention stability is provided after conditioning.

The combination of the robust YMC-Triart Diol-HILIC stationary phase with a cleaning step after each injection and an optimal equilibration step leads to a robust and reliable HILIC method.




Figure 1: Rapid conditioning of the YMC-Accura Triart Diol-HILIC column and the subsequent excellent reproducibility.






Figure 2 shows the non-targeted screening of a human plasma sample and selected chromatograms of polar metabolites with low concentrations of < 5 ppm.




Figure 2: Non-targeted HRMS analysis of human plasma: total ion current of the plasma sample (top) and selected ion chromatograms of annotated metabolites < 5 ppm (bottom).

The developed HILIC method using a bioinert YMC-Accura Triart Diol-HILIC column covers a wide range of polar compounds with an excellent peak capacity (baseline peak width of 0.25 min on average). In addition, the method achieves resolution of important critical pairs such as leucine and

isoleucine and asymmetric and symmetric dimethylarginine (ADMA and SDMA respectively). Together with the simultaneous high sensitivity, this method ensures a reliable generation of biological hypotheses.

* Application data by courtesy of Sergey Girel, Institute of Pharmaceutical Sciences of Western Switzerland (University of Geneva), Geneva, Switzerland.

