



Separation of plasmid isoforms using BioPro HIC BF

During bioprocessing of supercoiled DNA, chemical and physical factors can cause conformational changes leading to formation of the other, less favourable isoforms. Hydrophobic interaction chromatography (HIC) is a very

good method for monitoring the purity of DNA. In this application note, three plasmid DNA (pDNA) isoforms in linear, open circular and supercoiled conformation are separated by HIC using a BioPro HIC BF column.



The composition of the eluent is crucial for the successful separation of pDNA isoforms. Binding of pDNA to the stationary phase often depends on the correct amount of antichaotropic (kosmotropic) salt. In this application note this is achieved with 2.5M $(\text{NH}_4)_2\text{SO}_4$. To separate the three isoforms, two

buffers with different pH are tested. The separation of the pDNA isoforms is possible with both buffers (figure 1 & 2). By using ammonium phosphate buffer, the flow rate can be increased to 1.0mL/min, reducing the analysis time to approximately 19 min.



Table 1: Chromatographic conditions.

Column:	BioPro HIC BF (4µm) 100 x 4.6 mm ID
Part No.:	BHB00S04-1046WT
Eluents:	A) 50mM Tris-HCl containing 2.5 M (NH ₄) ₂ SO ₄ (pH 7.5) B) 50mM Tris HCl (pH 7.5)
Gradient:	4%B (0–3.5 min), 4–20%B (3.5–30.7 min), 20%B (30.7–37.7 min) 0.5
Flow rate:	mL/min
Temperature:	25°C
Injection:	9µL
Detection:	UV at 260nm
Sample:	Mixture of plasmid DNA (pUC19) isoforms in 2.0M (NH ₄) ₂ SO ₄ solution (open circular ¹ : 33µg/mL, linear ² : 33µg/mL, supercoiled ³ : 83µg/mL)

¹ pUC19 plasmid digested with nicking endonuclease Nt. BspQI

² pUC19 plasmid digested with restriction enzyme BamHI

³ pUC19 plasmid extracted from *E. coli* (2686 bp)




Figure 1: Separation of three pDNA isoforms open circular, linear and supercoiled using Tris-HCl (pH 7.5) as eluent.






Table 2: Chromatographic conditions.

Column:	BioPro HIC BF (4 μ m) 100 x 4.6 mm ID
Part No.:	BF00S04-1046WT
Eluents:	A) 50 mM NaH_2PO_4 - Na_2HPO_4 containing 2.5 M $(\text{NH}_4)_2\text{SO}_4$ (pH 6.5) B) 50 mM NaH_2PO_4 - Na_2HPO_4 (pH 6.5)
Gradient:	4%B (0–1.8 min), 4–20%B (1.8–15.4 min), 20%B (15.4–18.9 min)
Flow rate:	1.0 mL/min
Temperature:	40 °C
Injection:	9 μ L
Detection:	UV at 260 nm
Sample:	Mixture of plasmid DNA (pUC19) isoforms in 2.0 M $(\text{NH}_4)_2\text{SO}_4$ solution (open circular ¹ : 33 μ g/mL, linear ² : 33 μ g/mL, supercoiled ³ : 83 μ g/mL)

¹ pUC19 plasmid digested with nicking endonuclease Nt. BspQI

² pUC19 plasmid digested with restriction enzyme BamHI

³ pUC19 plasmid extracted from *E. coli* (2686 bp)




Figure 2: Separation of three pDNA isoforms open circular, linear and supercoiled using ammonium phosphate buffer (pH 6.5) as eluent.

