

NEW & UNIQUE

YMC CHIRAL ART

Immobilised Polysaccharides Phase


CHIRAL ART Cellulose-SJ



Features

- Unique Immobilised Polysaccharide Derivative
- Specifically Interesting for CHIRALPAK® OJ Users
- Stable in RP, NP and SFC Mode
- More Flexibility Due to Wide Range of Usable Solvents
- Highly Robust, also Suitable for SFC/SMB
- Extremely Attractive Pricing

Chiral Selector



Same chiral selector as in CHIRALPAK® OJ

Specifications


	CHIRAL ART Cellulose-SJ
Particle Size	3, 5 µm
Chiral Selector	Cellulose tris (4-methylbenzoate)
Type	Immobilised type
Separation Mode	Normal Phase / Reversed Phase / SFC
Usable pH range	2.0 – 9.0
Pressure limit	30 MPa (4350 psi)
(Alternative to)	Coated version CHIRALPAK® OJ)

now
available
3 µm

YMC
EUROPE GMBH

www.ymc.de

Unique Stereoselectivity Complements Other CHIRAL ART Phases



Ordering Information 5 μ m CHIRAL ART (Extract)

NEW

Dimension [mm]	Cellulose-SJ	Amylose-SA	Cellulose-SB	Cellulose-SC	Amylose-C	Cellulose-C
100 x 3.0	KSJ99S05-1003WT	KSA99S05-1003WT	KSB99S05-1003WT	KSC99S05-1003WT	KAN99S05-1003WT	KCN99S05-1003WT
100 x 4.6	KSJ99S05-1046WT	KSA99S05-1046WT	KSB99S05-1046WT	KSC99S05-1046WT	KAN99S05-1046WT	KCN99S05-1046WT
150 x 3.0	KSJ99S05-1503WT	KSA99S05-1503WT	KSB99S05-1503WT	KSC99S05-1503WT	KAN99S05-1503WT	KCN99S05-1503WT
150 x 4.6	KSJ99S05-1546WT	KSA99S05-1546WT	KSB99S05-1546WT	KSC99S05-1546WT	KAN99S05-1546WT	KCN99S05-1546WT
250 x 4.6	KSJ99S05-2546WT	KSA99S05-2546WT	KSB99S05-2546WT	KSC99S05-2546WT	KAN99S05-2546WT	KCN99S05-2546WT

Several other dimensions or particle sizes available. Check www.ymc.de/chiral-columns.html

YMC CO., LTD.
YMC Karasuma-Gojo Bld. 284 Daigo-cho,
Karasuma Nishiru Gojo-dori Shimogyo-ku,
Kyoto 600-8106 Japan
Phone +81(0)75-342-4515, FAX +81(0)75-342-4550
www.ymc.co.jp

YMC Co., Ltd. Shanghai Rep. Office
www.ymchina.com

YMC Europe GmbH
Schöttmannshof 19
D-46539 Dinslaken
Germany
Phone +49(0)2064/427-0, FAX +49(0)2064/427-222
www.ymc.de

YMC India Pvt. Ltd.
www.ymcindia.com

YMC Korea Co., Ltd.
www.ymckorea.com

YMC Schweiz GmbH
Im Wasenboden 8
4056 Basel
Phone +41 61 561 80 50, Fax +41 61 561 80 59
www.ymc-schweiz.ch

YMC Taiwan Co., Ltd.
www.ymctaiwan.com

YMC America, Inc.
941 Marcon Boulevard Suite 201
Allentown, PA18109 USA
Phone +1-610-266-8650, FAX +1-610-266-8652
www.ymcamerica.com

YMC Singapore Tradelinks Pte. Ltd.
www.ymc.sg