



\* Application Data by courtesy YMC America, Inc.

Column: Meteoric Core C18, 2.7  $\mu\text{m}$ , 100 x 4.6 mm ID  
Part No.: CAS08SQ7-1046PT  
Eluent: water+formic acid (0.1%) / acetonitrile (27/73)  
Flow rate: 1.25 ml/min  
Detection: UV at 220 nm  
Injection: 25  $\mu\text{L}$  (0.2 mg/mL)  
Temperature: 35  $^{\circ}\text{C}$