



Column : YMC-Pack Diol-200 + Diol-120  
 500 × 8.0 mm I.D. × 2  
 Eluent : 0.1M KH<sub>2</sub>PO<sub>4</sub> - K<sub>2</sub>HPO<sub>4</sub> (pH 7.0) containing 0.2M NaCl  
 Flow rate : 1.0 mL/min  
 Temperature : ambient (24 °C)  
 Detection : UV at 215 nm, 0.08 AUFS  
 Injection : 30 μL (0.2 mg/mL)