



1. Fast Green FCF  
(FD&C Green NO.3, C.I.42053, G-3)
2. Brilliant Blue FCF  
(FD&C Blue No.1, C.I.42090, B-1)
3. Erythrosine  
(FD&C Red No.3, C.I.45430, R-3)
4. Acid Red  
(C.I.45100, R-106)
5. Phloxine  
(C.I.45410, R-104)
6. Rose Bengale  
(C.I.45440, R-105)

Column : YMC-Pack ODS-A (5  $\mu$ m, 12nm)  
 150  $\times$  4.6 mm I.D.  
 Eluent : acetonitrile/10mM CH<sub>3</sub>COONH<sub>4</sub> (25/75)  
 Flow rate : 1.0 mL/min  
 Temperature : 30  $^{\circ}$ C  
 Detection : UV at 254 nm, 0.08 AUFS  
 Injection : 10  $\mu$ L (35~140  $\mu$ g/mL)