



- 1. 5'-CMP
- 2. 5'-AMP
- 3. 5'-UMP
- 4. 5'-IMP
- 5. 5'-GMP

Column : YMC-Pack Polyamine II
 250 × 4.6 mm I.D.
 Eluent : 50mM KH₂PO₄ - H₃PO₄ (pH 3.5)
 Flow rate : 1.0 mL/min
 Temperature : 40 °C
 Detection : UV at 260 nm, 0.16 AUFS
 Injection : 10 μL (0.05~0.1 mg/mL)