



1. Acetaminophene
2. 2-Acetamidophenol
3. Caffeine
4. Acetanilide
5. Acetylsalicylic Acid
6. Phenacetine
7. Salicylic Acid

Column: YMC-Triart C18, 1.9  $\mu\text{m}$ , 12 nm (50 x 2.0 mm ID)  
Part-No.: TA12SP9-0502PT  
Eluent: water + formic acid (pH 2.5) / acetonitrile (50/50)  
Flow rate: 0.8 ml/min  
Detection: UV at 240 nm  
Injection: 1  $\mu\text{l}$   
Temperature: 40  $^{\circ}\text{C}$