



1. Resorcinol
2. Pyrocatechol
3. 2-Methylresorcinol
4. 4-Methylcatechol
5. 2,5-Methylresorcinol
6. 4-Nitrocatechol

Column: YMC-Triart PFP, 1.9 μm (100 x 2.0 mm ID)
Part-No.: TPF12SP9-1002PT
Eluent: water + 0.1% formic acid / acetonitrile + 0.1% formic acid (85/15)
Flow rate: 0.8 ml/min
Detection: UV at 270 nm
Injection: 0.5 μl
Temperature: 25 $^{\circ}\text{C}$