

# YMC Recommended Phase Alternatives by Manufacturer



Column	YMC Alternative	YMC Recommended Alternative
<b>Agilent</b>		
AdvanceBio HIC	–	BioPro HIC HT
AdvanceBio RP-mAb SB-C8	–	YMC-Triart Bio C4
AdvanceBio RP-mAb C4	–	YMC-Triart Bio C4
AdvanceBio RP-mAb diphenyl	–	YMC-Triart Bio C4
AdvanceBio Amino Acid Analysis	–	Meteoric Core C18
AdvanceBio SEC - 130	–	YMC-Pack Diol-120
AdvanceBio SEC - 200	YMC-Pack Diol-200	YMC-SEC MAB
AdvanceBio SEC - 300	YMC-Pack Diol-300	
ProSEC 300S	YMC-Pack Diol-300	
Bio-SCX	BioPro IEX SF	
Bio-SAX	BioPro IEX QF	
Bio-WCX	–	BioPro IEX SF
Bio-WAX	–	BioPro IEX QF
Bio Mab NP	–	BioPro IEX SF
Polaris C18-A	YMC-Pack ODS-A	YMC-Triart C18
Polaris C8-A	–	YMC-Triart C8
Polaris C18-Ether	–	YMC-Triart C18
Polaris C8-Ether	–	YMC-Triart C8
Polaris Amide-C18	–	YMC-Triart C18
Polaris NH2	YMC-Pack NH2	
Polaris Si-A	YMC-Pack SIL	
Poroshell 120 EC-C18	Meteoric Core C18	YMC-Triart C18
Poroshell 120 EC-C8	Meteoric Core C8	YMC-Triart C8
Poroshell 120 HILIC	–	YMC-Triart Diol-HILIC
Poroshell 120 SB-C18	–	YMC-Triart C18
Poroshell 120 SB-C8	–	YMC-Triart C8
Poroshell 120 Phenyl-Hexyl	–	YMC-Triart Phenyl
Poroshell 120 Bonus-RP	–	Meteoric Core C18
Poroshell 120 PFP	–	YMC-Triart PFP
Poroshell HPH-C18	–	YMC-Triart C18
Poroshell HPH-C8	–	YMC-Triart C8
Poroshell 120 EC-CN	–	YMC-Pack CN
Poroshell 300 SB-C3	–	YMC-Triart Bio C4
Poroshell 300 SB-C8	–	YMC-Pack C8
Poroshell 300 SB-C18	–	YMC-Triart Bio C18
Poroshell 300 Extend-C18	–	YMC-Triart Bio C18
Pursuit C18	YMC-Pack ODS-A	YMC-Triart C18
Pursuit C8	YMC-Pack Octyl	YMC-Triart C8
Pursuit PAH	YMC PAH	
Pursuit Diphenyl	–	YMC-Triart Phenyl
Pursuit PFP	–	YMC-Triart PFP
ZORBAX Eclipse Plus C18	YMC-Pack Pro C18	YMC-Triart C18
ZORBAX Eclipse Plus C8	YMC-Pack Pro C8	YMC-Triart C8
ZORBAX Eclipse Phenyl-Hexyl	-	YMC-Triart Phenyl
ZORBAX Eclipse PAH	YMC PAH	
ZORBAX SB-C18	–	YMC-Pack ODS-AL
ZORBAX SB-C8	–	YMC-Triart C8
ZORBAX SB-Phenyl	–	YMC-Triart Phenyl
ZORBAX SB-CN	–	YMC-Pack CN
ZORBAX SB-Aq	–	YMC-Pack ODS-AQ
ZORBAX Eclipse XDB-C18	YMC-Pack Pro C18	YMC-Triart C18

Column	YMC Alternative	YMC Recommended Alternative
ZORBAX Eclipse XDB-C8	YMC-Pack Pro C8	YMC-Triart C8
ZORBAX Eclipse XDB-Phenyl	YMC-Pack Phenyl	YMC-Triart Phenyl
ZORBAX Eclipse XDB-CN	YMC-Pack CN	
ZORBAX Extend C18	–	YMC-Triart C18
ZORBAX Bonus-RP	–	YMC-Triart C18
ZORBAX Rx-C18	–	YMC-Pack ODS-AL
ZORBAX Rx-C8	–	YMC-Triart C8
ZORBAX HILIC Plus	–	YMC-Triart Diol-HILIC
ZORBAX NH2	YMC-Pack NH2	
Zorbax 300Å Extend-C18	–	YMC-Triart Bio C18
ZORBAX GF-250	–	YMC-Pack Diol-200
ZORBAX GF-450	YMC-Pack Diol-300	
<b>Daicel</b>		
Chiralpak AD(-3/-H)	CHIRAL ART Amylose-C	CHIRAL ART Amylose-C Neo
Chiralcel OD(-3/-H)	CHIRAL ART Cellulose-C	
Chiralpak OJ(-3/-H)	–	CHIRAL ART Cellulose-SJ
Chiralpak IA(-3/-U)	CHIRAL ART Amylose-SA	
Chiralpak IB(-3/-U)	CHIRAL ART Cellulose-SB	
Chiralpak IC(-3/-U)	CHIRAL ART Cellulose-SC	
<b>Phenomenex</b>		
Aeris PEPTIDE XB-C18	YMC-Triart C18	
Aeris WIDEPOR C4	YMC-Triart Bio C4	
Aeris WIDEPOR XB-C18	YMC-Triart Bio C18	
Aeris WIDEPOR XB-C8	YMC-Pack C8, 200 Å	YMC-Triart Bio C4
bioZen Peptide PS-C18	YMC-Triart C18	
bioZen Peptide XB-C18	YMC-Triart C18	
bioZen Intact C4	YMC-Triart Bio C4	
bioZen Intact XB-C8	–	YMC-Triart Bio C4
bioZen WCX	–	BioPro IEX SF
bioZen SEC-2	–	YMC-Pack Diol-200
bioZen SEC-3	YMC-Pack Diol-300	YMC-SEC MAB
Clarity Oligo-SAX	–	BioPro IEX QF
Clarity Oligo-MS	YMC-Triart C18	
Clarity Oligo-RP	YMC-Triart C18	
Clarity Oligo-XT	YMC-Triart C18	
Clarity Oligo-WAX	–	BioPro IEX QF
Gemini C18	YMC-Triart C18	
Gemini C6-Phenyl	–	YMC-Triart Phenyl
Gemini NX-C18	YMC-Triart C18	
Jupiter C4	YMC-Pack C4, 300 Å	YMC-Triart Bio C4
Jupiter C18	YMC-Pack ODS-A, 300 Å	YMC-Triart Bio C18
Kinetex EVO C18	YMC-Triart C18	
Kinetex XB C18	YMC-Triart C18	
Kinetex PS C18	YMC-Triart C18	
Kinetex Polar C18	YMC-Triart C18	
Kinetex C18	YMC-Triart C18	
Kinetex Phenyl-Hexyl	YMC-Triart Phenyl	
Kinetex HILIC	YMC-Triart Diol-HILIC	
Kinetex PFP	YMC-Triart PFP	
Kinetex C8	YMC-Triart C8	
Kinetex F5	YMC-Triart PFP	
Kinetex PAH	YMC-PAH	

# YMC Recommended Phase Alternatives by Manufacturer



Column	YMC Alternative	YMC Recommended Alternative
Luna C18	YMC-Pack Pro C18	YMC-Triart C18
Luna C18(2)	YMC-Pack Pro C18	YMC-Triart C18
Luna C8	YMC-Pack Pro C8	YMC-Triart C8
Luna C8(2)	YMC-Pack Pro C8	YMC-Triart C8
Luna CN	YMC-Pack CN	
Luna HILIC	YMC-Pack Diol-NP	YMC-Triart Diol-HILIC
Luna NH2	YMC-Pack NH2	
Luna Omega C18	YMC-Triart C18	
Luna Omega Polar C18	Hydrosphere C18 (for HPLC)	
Luna Omega PS C18	–	YMC-Triart C18
Luna PFP (2)	–	YMC-Triart PFP
Luna Phenyl-Hexyl	–	YMC-Triart Phenyl
Luna SCX	–	BioPro IEX SF
Luna Silica(2)	YMC-Pack SIL	
Lux Amylose-1	CHIRAL ART Amylose-C	
Lux i-Amylose-1	CHIRAL ART Amylose-SA	
Lux Cellulose-1	CHIRAL ART Cellulose-C	
Lux i-Cellulose-5	CHIRAL ART Cellulose-SC	
Lux Cellulose-3	–	CHIRAL ART Cellulose-SJ
Synergi Fusion-RP	–	YMC-Triart C18
Synergi Hydro-RP	Hydrosphere C18	YMC-Triart C18
Yarra SEC-2000	–	YMC-Pack Diol-200
Yarra SEC-3000	YMC-Pack Diol-300	YMC-SEC MAB
Yarra SEC-4000	–	YMC-Pack Diol-300
Yarra SEC-X150	–	YMC-Pack Diol-200
Yarra SEC-X300	YMC-Pack Diol-300	YMC-SEC MAB
<b>Sigma-Adrich + Merck Millipore</b>		
apHera Amino	–	YMC-Pack Polyamine II
apHera C18	–	YMC-Pack Polymer C18
apHera C8	–	YMC-Pack C8 300 Å
apHera C4	–	YMC-Triart Bio C4
Ascentis C18	YMC-Pack Pro C18	YMC-Triart C18
Ascentis C8	YMC-Pack Pro C8	YMC-Triart C8
Ascentis ES-Cyano	YMC-Pack CN	
Ascentis Silica	YMC-Pack SIL	
Ascentis Phenyl	YMC-Triart Phenyl	
Ascentis Express C18	Meteoric Core C18	YMC-Triart C18
Ascentis Express AQ-C18	–	YMC-Triart C18
Ascentis Express C8	Meteoric Core C8	YMC-Triart C8
Ascentis Express C30	–	YMC Carotenoid
Ascentis Express ES-Cyano	–	YMC-Pack CN
Ascentis Express F5 (PFP)	–	YMC-Triart PFP
Ascentis Express Phenyl-Hexyl	–	YMC-Triart Phenyl
Ascentis Express HILIC (Si)	–	YMC-Triart Diol-HILIC
Astec CYCLOBOND I 2000	–	YMC CHIRAL beta-CD BR
Astec CYCLOBOND II	–	YMC CHIRAL gamma-CD BR
Astec Cellulose DMP	CHIRAL ART Cellulose-C	
BIOshell 400 Å Protein C4	–	YMC-Triart Bio C4
BIOshell 400 Å Protein C18	–	YMC-Triart Bio C18
BIOshell 160 Å Peptide C18	Meteoric Core C18 Bio	YMC-Triart C18
BIOshell 160 Å Peptide CN	–	YMC-Pack CN
BIOshell 160 Å Peptide Phenyl-Hexyl	–	YMC-Triart Phenyl

Column	YMC Alternative	YMC Recommended Alternative
Discorvery BIOWide Pore C18	YMC-Pack ODS-A 300 Å	YMC-Triart Bio C18
Discorvery BIOWide Pore C8	YMC-Pack C8 300 Å	
Discorvery BIOWide Pore C5	–	YMC-Triart Bio C4
Discovery C18	YMC-Pack ODS-A	YMC-Triart C18
Discovery C8	YMC-Pack C8	YMC-Triart C8
Discovery Cyano	YMC-Pack Diol-NP	
Discovery HS C18	YMC-Pack ODS-A	YMC-Triart C18
Discovery HS F5	–	YMC-Triart PFP
LiChrosorb RP-18	YMC-Pack Pro C18	YMC-Triart C18
LiChrosorb RP-8	YMC-Pack Pro C8	YMC-Triart C8
LiChrosorb Si 100	YMC-Pack SIL	
LiChrospher 100 CN	YMC-Pack CN	
LiChrospher 100 Diol	YMC-Pack Diol	YMC-Triart Diol-HILIC
LiChrospher 100 NH2	YMC-Pack NH2	
LiChrospher 100 RP-18	YMC-Pack ODS-AL	YMC-Triart C18
LiChrospher 100 RP-18 endcapped	YMC-Pack Pro C18	YMC-Triart C18
LiChrospher 100 RP-8	–	YMC-Pack Pro C8
LiChrospher 100 RP-8 endcapped	YMC-Pack Pro C8	YMC-Triart C8
LiChrospher Si 100	YMC-Pack SIL	
LiChrospher Si 60	YMC-Pack SIL	
Purospher RP-18	YMC-Pack Pro C18 RS	YMC-Triart C18 ExRS
Purospher RP-18 endcapped	YMC-Pack Pro C18 RS	YMC-Triart C18 ExRS
Purospher RP-18 HC	–	YMC-Pack ODS-AL
Purospher STAR NH2	YMC-Pack NH2	
Purospher STAR RP-18 endcapped	YMC-Pack ODS-AQ	YMC-Triart C18
Purospher STAR RP-18 endcapped UHPLC	–	YMC-Triart C18
Purospher STAR RP-8 endcapped	YMC-Pack Pro C8	YMC-Triart C8
Purospher STAR Si	YMC-Pack SIL	
Superspher 100 RP-18	YMC-Pack Pro C18	YMC-Triart C18
Superspher 100 RP-18 endcapped	YMC-Pack Pro C18	YMC-Triart C18
Superspher Si 60	YMC-Pack SIL	
SUPELCOSIL LC-C18	YMC-Pack ODS-A	
SUPELCOSIL LC-C18-DB	YMC-Pack ODS-A	
SUPELCOSIL LC-C8	YMC-Pack C8	
SUPELCOSIL LC-C8-DB	YMC-Pack C8	
SUPELCOSIL LC-1	YMC-Pack TMS	
SUPELCOSIL LC-CN	YMC-Pack CN	
SUPELCOSIL LC-PCN	YMC-Pack CN	
SUPELCOSIL LC-Diol	YMC-Pack Diol-NP	
SUPELCOSIL LC-NH2	YMC-Pack NH2	
SUPELCOSIL LC-Si	YMC-Pack Sil	
Titan C18 UHPLC	YMC-Triart C18	
<b>Thermo Fisher</b>		
Acclaim 120 C18	YMC-Pack Pro C18	YMC-Triart C18
Acclaim HILIC-1	–	YMC-Triart Diol-HILIC
Acclaim PepMap 100-C18	YMC-Pack Pro C18	
Acclaim VANQUISH	YMC-Pack Pro C18	
Acclaim C30	YMC Carotenoid	
Acclaim Phenyl-1	YMC-Pack Phenyl	YMC-Triart Phenyl
Acclaim HILIC-10	–	YMC-Triart Diol-HILIC
Acclaim 120 C8	YMC-Pack Pro C8	YMC-Triart C8
Acclaim SEC-300	YMC-Pack Diol-300	

# YMC Recommended Phase Alternatives by Manufacturer



Column	YMC Alternative	YMC Recommended Alternative
Acclaim 300-C18	YMC-Pack ODS-A	YMC-Triart Bio C18
Accucore 150 C18	Meteoric Core C18 BIO	
Accucore C18	Meteoric Core C18	YMC-Triart C18
Accucore C30	–	YMC Carotenoid
Accucore C8	Meteoric Core C8	YMC-Triart C8
Accucore 150 C4	–	YMC-Pack Pro C4
Accucore HILIC	–	YMC-Triart Diol-HILIC
Accucore PFP	–	YMC-Triart PFP
Accucore Phenyl-Hexyl	–	YMC-Triart Phenyl
Accucore Polar Premium	–	Meteoric Core C18 BIO
Accucore RP-MS	Meteoric Core C18	YMC-Triart C18
Accucore Vanquish C18+ UHPLC	–	YMC-Triart C18
Accucore aQ C18 Polar Endcapped	–	YMC-Triart C18
Accucore Phenyl-X	–	YMC-Pack Phenyl
BioBasic 18	YMC-Pack ODS-A	YMC-Triart Bio C18
BioBasic 8	YMC-Pack C8	YMC-Triart Bio C4
BioBasic 4	YMC-Pack C4	YMC-Triart Bio C4
BioBasic AX	–	BioPro IEX QA
BioBasic SCX	BioPro IEX SP	
BioBasic SEC 60	YMC-Pack Diol-60	
BioBasic SEC 120	YMC-Pack Diol-120	
BioBasic SEC 300	YMC-Pack Diol-300	YMC-SEC MAB
Hypersil Chiral AT	CHIRAL ART Amylose-C	
Hypersil Chiral OT	CHIRAL ART Cellulose-C	
Hypersil GOLD Amino	–	YMC-Pack NH2
Hypersil GOLD aQ	YMC-Pack ODS-AQ	YMC-Triart C18
Hypersil GOLD C18	YMC-Pack Pro C18	YMC-Triart C18
Hypersil GOLD C4	YMC-Pack C4	YMC-Pack Pro C4
Hypersil GOLD C8	YMC-Pack Pro C8	YMC-Triart C8
Hypersil GOLD Cyano	YMC-Pack CN	
Hypersil GOLD PFP	–	YMC-Triart PFP
Hypersil GOLD Phenyl	–	YMC-Triart Phenyl
Hypersil GOLD SAX	–	BioPro IEX QA
Hypersil GOLD AX	–	BioPro IEX QA
Hypersil GOLD Silica	YMC-Pack SIL	
Hypersil GOLD PEI HILIC	–	YMC-Triart Diol-HILIC
Hypersil BDS C18	YMC-Pack Pro C18	YMC-Triart C18
Hypersil BDS C8	YMC-Pack Pro C8	YMC-Triart C8
Hypersil BDS Phenyl	YMC-Pack Phenyl	YMC-Triart Phenyl
Hypersil BDS Cyano	YMC-Pack CN	
Hypersil ODS	YMC-Pack ODS-A	YMC-Triart C18
Hypersil ODS-2	YMC-Pack ODS-A	YMC-Triart C18
Hypersil Mos-2	YMC-Pack C8	YMC-Triart C8
Hypersil SAS	YMC-Pack TMS	
Hypersil Phenyl-2	YMC-Pack Phenyl	YMC-Triart Phenyl
Hypersil CPS-2	YMC-Pack CN	
Hypersil APS-2	YMC-Pack NH2	
Hypersil Silica	YMC-Pack SIL	
Hypersil Green PAH	YMC-Pack PAH	
MABPac SCX-10	BioPro IEX SF	
MabPac SEC-1	YMC-Pack Diol-300	
MABPac HIC-Butyl	BioPro HIC HT	

Column	YMC Alternative	YMC Recommended Alternative
MABPac HIC-10	–	BioPro HIC HT
MABPac HIC-20	–	BioPro HIC HT
ProPac WCX-10	–	BioPro IEX SF
ProPac SAX-10	–	BioPro IEX QF
ProPac SCX-10	–	BioPro IEX SF
ProPac PA1	–	BioPro IEX QF
ProPac SCX-20	–	BioPro IEX SF
ProPac WAX-10	–	BioPro IEX QF
ProPac Elite	–	BioPro IEX SF
ProSwift RP-4H	–	YMC-Triart Phenyl
Synchromis Amino	YMC-Pack NH2	
Synchromis aQ C18	YMC-Pack ODS-AQ	YMC-Triart C18
Synchromis C18	YMC-Pack Pro C18	YMC-Triart C18
Synchromis C8	YMC-Pack Pro C8	YMC-Triart C8
Synchromis Phenyl	YMC-Pack Phenyl	YMC-Triart Phenyl
Synchromis Silica	YMC-Pack SIL	YMC-Pack SIL
Synchromis HILIC	–	YMC-Triart Diol-HILIC
<b>Tosoh</b>		
TSKgel CN-80Ts	–	YMC-Pack CN
TSKgel Octyl-80Ts	YMC-Pack C8	YMC-Triart C8
TSKgel ODS-80Tm	YMC-Pack Pro C18	YMC-Triart C18
TSKgel ODS-80Ts	YMC-Pack Pro C18	YMC-Triart C18
TSKgel ODS-80Ts QA	YMC-Pack Pro C18	YMC-Triart C18
TSKgel ODS-100V	YMC-Pack ODS-AQ	YMC-Triart C18
TSKgel ODS-100Z	YMC-Pack Pro C18	YMC-Triart C18
TSKgel Super-Octyl	Meteoric Core C8	YMC-Triart C8
TSKgel Super-ODS	Meteoric Core C18	YMC-Triart C18
TSKgel Super-Phenyl	–	YMC-Triart Phenyl
TSKgel ODS-120A	–	YMC-Triart C18
TSKgel ODS-120T	–	YMC-Triart C18
TSKgel OligoDNA-RP	–	YMC-Triart Bio C18
TSKgel ODS-140HTP	–	Hydrosphere C18
TSKgel Protein C4-300	YMC-Pack C4, 300 Å	YMC-Triart Bio C4
TSKgel TMS-250	–	YMC-Pack TMS
TSKgel NH2-100	YMC-Pack NH2	
TSKgel DNA-STAT	BioPro IEX QF	
TSKgel BioAssist Q	BioPro IEX QA	
TSKgel BioAssist S	BioPro IEX SP	
TSKgel Butyl-NPR	BioPro HIC HT	
TSKgel SP-NPR	–	BioPro IEX SF
TSKgel SuperSW mAb HR	–	YMC-Pack Diol-300
TSKgel SuperSW mAb HTP	–	YMC-Pack Diol-300
TSKgel Ultra SW Aggregate	YMC-Pack Diol-300	
TSKgel G2000SW	–	YMC-Pack Diol-120
TSKgel G3000SW	–	YMC-Pack Diol-300
TSKgel G2000WXL	–	YMC-Pack Diol-120
TSKgel G3000SWXL	–	YMC-Pack Diol-300
TSKgel SuperSW2000	–	YMC-Pack Diol-120
TSKgel SuperSW3000	YMC-SEC MAB	
TSKgel UP-SW3000	–	YMC-SEC MAB

# YMC Recommended Phase Alternatives by Manufacturer



Column	YMC Alternative	YMC Recommended Alternative
<b>Waters</b>		
ACQUITY UPLC BEH C18	YMC-Triart C18	
ACQUITY UPLC BEH C8	YMC-Triart C8	
ACQUITY UPLC BEH Phenyl	YMC-Triart Phenyl	
ACQUITY UPLC BEH Shield-RP18	–	YMC-Triart C18
ACQUITY UPLC BEH HILIC	–	YMC-Triart Diol-HILIC
ACQUITY UPLC BEH Amide	–	YMC-Triart Diol-HILIC
ACQUITY UPLC BEH300 C18	YMC-Triart Bio C18	
ACQUITY UPLC BEH300 C4	YMC-Triart Bio C4	
ACQUITY UPLC OST C18	YMC-Triart C18	
ACQUITY UPLC BEH450 SEC	–	YMC-Pack Diol-300
ACQUITY UPLC BEH200 SEC	–	YMC-SEC MAB
ACQUITY UPLC BEH125 SEC	–	YMC-Pack Diol-120
ACQUITY UPLC CSH C18	–	YMC-Triart C18
ACQUITY UPLC CSH Phenyl-Hexyl	–	YMC-Triart Phenyl
ACQUITY UPLC CSH Fluoro-Phenyl	–	YMC-Triart PFP
ACQUITY UPLC HSS-T3	–	YMC-Triart C18
ACQUITY UPLC HSS-C18	–	YMC-Triart C18
ACQUITY UPLC HSS-PFP	–	YMC-Triart PFP
ACQUITY UPLC BEH Peptide C18	YMC-Triart C18	
ACQUITY UPLC Peptide BEH C4	YMC-Triart Bio C4	
ACQUITY UPLC Peptide CSH C18	–	YMC-Triart C18
ACQUITY UPLC Peptide HSS T3	YMC-Triart C18	
ACQUITY UPLC Protein BEH C4	YMC-Triart Bio C4	
ACQUITY UPLC Protein BEH SEC 125	–	YMC-Pack Diol-120
ACQUITY UPLC Protein BEH SEC 200	–	YMC-SEC MAB
ACQUITY UPLC Protein BEH SEC 450	–	YMC-Pack Diol-300
Atlantis dC18	–	YMC-Triart C18
Atlantis HILIC	YMC-Pack SIL	YMC-Triart Diol-HILIC
Atlantis T3	YMC-Pack ODS-AQ	YMC-Triart C18
Atlantis PREMIER BEH C18 AX	–	YMC-Triart C18
XBridge BEH C18	YMC-Triart C18	
XBridge BEH C8	YMC-Triart C8	
XBridge BEH Shield RP18 XP	YMC-Triart C18	
XBridge BEH Phenyl	YMC-Triart Phenyl	
XBridge BEH HILIC	–	YMC-Triart Diol-HILIC
XBridge BEH Amide	–	YMC-Triart Diol-HILIC
XBridge Peptide BEH C18 130A	YMC-Triart C18	
XBridge Peptide BEH C18 300A	YMC-Triart Bio C18	
XBridge Protein BEH C4 300A	YMC-Triart Bio C4	
XBridge Oligonucleotide BEH C18	YMC-Triart C18	
XBridge Protein BEH SEC 125	–	YMC-Pack Diol-120
XBridge Protein BEH SEC 200	–	YMC-SEC MAB
XBridge Protein BEH SEC 450	–	YMC-Pack Diol-300
BioResolve RP mAb	–	YMC-Triart Bio C4
BioResolve SCX mAb	BioPro IEX SF	
BioSuite CM	–	BioPro IEX SP
BioSuite Q-Peek	–	BioPro IEX QA
BioSuite SP	–	BioPro IEX SP
BioSuite SP-Peek	–	BioPro IEX SP
BioSuite C18	–	YMC-Triart C18

Column	YMC Alternative	YMC Recommended Alternative
BioSuite C18 300	–	YMC-Triart Bio C18
BioSuite DEAE	–	BioPro IEX QA
BioSuite Q	–	BioPro IEX QA
BioSuite HR SEC 125	YMC-Pack Diol-120	
BioSuite HR SEC 250	YMC-SEC MAB	
BioSuite HR SEC 450	–	YMC-Pack Diol-300
BioSuite SEC 125	–	YMC-Pack Diol-120
BioSuite SEC 250	–	YMC-SEC MAB
BioSuite SEC 450	–	YMC-Pack Diol-300
BioSuite UHR SEC 125	–	YMC-Pack Diol-120
BioSuite UHR SEC 250	–	YMC-SEC MAB
Cortecs C18+	–	Meteoric Core C18
Cortecs C18	Meteoric Core C18	YMC-Triart C18
Cortecs HILIC	–	YMC-Triart Diol-HILIC
Cortecs C8	Meteoric Core C8	YMC-Triart C8
Cortecs Phenyl	–	YMC-Triart Phenyl
Delta-Pak C18	YMC-Triart C18	
Delta-Pak C4, 100 A	YMC-Pack C4	
Delta-Pak C4, 300 A	YMC-Triart Bio C4	
μBondapak	–	YMC-Triart C18
PROTEIN-PAK DEAE	–	BioPro IEX QA
Protein-Pak Hi Res CM	BioPro IEX QF	
Protein-Pak Hi Res HIC	BioPro HIC HT	
Protein-Pak Hi Res Q	BioPro IEX QF	
Protein-Pak Hi Res SP	BioPro IEX QF	
PROTEIN-PAK SEC 60	–	YMC-Pack Diol-60
PROTEIN-PAK SEC 125	–	YMC-Pack Diol-120
PROTEIN-PAK SEC 300	–	YMC-SEC MAB
SunFire C18	YMC-Pack Pro C18	YMC-Triart C18
SunFire C8	YMC-Pack Pro C8	YMC-Triart C8
SunFire Silica	YMC-Pack SIL	
Symmetry C18	YMC-Pack Pro C18	YMC-Triart C18
Symmetry C8	YMC-Pack Pro C8	YMC-Triart C8
Symmetry300 C18	YMC-Pack ODS-A, 300 Å	YMC-Triart Bio C18
Symmetry300 C4	YMC-Pack C4, 300 Å	YMC-Triart Bio C4
SymmetryShield RP18	–	YMC-Triart C18
SymmetryShield RP8	–	YMC-Triart C8
SymmetryPrep C18	–	Triart Prep C18
SymmetryPrep C8	–	Triart Prep C8
Torus Diol	–	YMC-Triart Diol
Viridis CSH Fluoro-Phenyl	–	YMC-Triart PFP
Viridis HSS C18 SB	–	YMC-Triart C18
Viridis Silica	–	YMC-Pack SIL
XSelect Peptide BEH C18	YMC-Triart C18	
XSelect Peptide CSH C18	YMC-Triart C18	
XSelect Peptide HSS T3	YMC-Triart C18	
XTerra MS C18	–	YMC-Triart C18
XTerra MS C8	–	YMC-Triart C8
XTerra Phenyl	–	YMC-Triart Phenyl
XTerra RP18	–	YMC-Triart C18
XTerra RP8	–	YMC-Triart C8