



# Easy ProFamily Column Selection

for Method Development




## YMC-Pack Pro C18




pH-stable 2–8  
ideal for basic, acidic and polar compounds

## YMC-Pack Pro C8




pH-stable 2.0–7.5  
good alternative to C18 phases

## YMC-Pack Pro C18 RS




pH-stable 1–10  
hydrophobic compounds  
steric selectivity

## YMC-Pack Pro C4



pH-stable 2.0–7.5  
ideal for faster separations

## Hydrosphere C18



pH-stable 2–8  
also 100% aqueous eluents



Schöttmannshof 19 · D-46539 Dinslaken  
Phone: +49 2064 427-0 · Fax: +49 2064 427-222 · Email info@ymc.eu

[www.ymc.eu](http://www.ymc.eu)

