



YMC-Triart C18




versatile applications
first choice for
method development
pH-stable 1–12
also 100% aqueous
eluents

YMC-Triart C18 ExRS



pH-stable 1–12
hydrophobic substances
positional isomers


YMC-Triart C8




alternative to C18
short retention time
pH-stable 1–12

Easy YMC-Triart Column Selection

for Method Development




YMC-Triart Phenyl




aromatic compounds
(π -electron acceptor)
conjugated systems
also 100% aqueous
eluents

YMC-Triart PFP



aromatic compounds
(π -electron donor)
cis-trans isomers
polar halogenated
compounds
also 100% aqueous
eluents



YMC
EUROPE GMBH

Schöttmannshof 19 · D-46539 Dinslaken
Phone: +49 2064 427-0 · Fax: +49 2064 427-222 · Email info@ymc.eu

www.ymc.eu