

**NEW &
INNOVATIVE**

YMC CHIRAL ART

Immobilised Polysaccharide Phase
CHIRAL ART Cellulose-SZ



Benefit from superior quality and service for your chiral separation needs!

- Free chiral screening service for your compounds
- Expert application and troubleshooting support
- Contract method development and purification service
- Screening kits of 3 different CHIRAL ART columns available at special price

Enhance your chiral portfolio with robust and reproducible columns:

Product Name	Particle Size	CHIRAL Selector	Type	Competitive Products	
				CHIRALPAK®/CHIRALCEL®	Lux®
CHIRAL ART Amylose-C (Neo)	3 µm 5 µm 10 µm 20 µm	Amylose tris (3,5-dimethylphenylcarbamate)	coated	AD(-H/-3)	Amylose-1
CHIRAL ART Cellulose-C		Cellulose tris (3,5-dimethylphenylcarbamate)		OD(-H/-3)	Cellulose-1
CHIRAL ART Amylose-SA		Amylose tris (3,5-dimethylphenylcarbamate)	immobilised	IA(-3)	i-Amylose-1
CHIRAL ART Cellulose-SB		Cellulose tris (3,5-dimethylphenylcarbamate)		IB(-3)	n.a.
CHIRAL ART Cellulose-SC		Cellulose tris (3,5-dichlorophenylcarbamate)		IC(-3)	i-Cellulose-5
CHIRAL ART Cellulose-SJ		Cellulose tris (4-methylbenzoate)		IJ(-3)	[coated Cellulose-3]
NEW CHIRAL ART Cellulose-SZ	3 µm 5 µm	Cellulose tris (3-chloro-4-methylphenylcarbamate)		[coated OZ]	[coated Cellulose-2]

YMC CO., LTD.

YMC Karasuma-Gojo Bld. 284 Daigo-cho,
Karasuma Nishiiru Gojo-dori Shimogyo-ku,
Kyoto 600-8106 Japan
Phone +81(0)75-342-4515, FAX +81(0)75-342-4550
www.ymc.co.jp

YMC Europe GmbH

Schöttmannshof 19
46539 Dinslaken
Germany
Phone +49(0)2064/427-0, FAX +49(0)2064/427-222
www.ymc.de

YMC Schweiz GmbH

Im Wasenboden 8
4056 Basel
Phone + 41 61 561 80 50, Fax + 41 61 561 80 59
www.ymc-schweiz.ch