


Accessories for YMC Glass Columns




Accessories-TAC05

Accessories for ECO^{PLUS} glass columns with 5 mm ID.

Connection Principle (Tubing is pre-assembled):




Connection to
System
1/4"-28G



JR-CFL-CB1PK-S

JR-55050

Connection to
System
10-32 HPLC



JR-ZC1FPKB

JR-58000

Connection to
System
M6




K.P200N

K.P207

Accessories-1/16

Accessories for ECO and ECO^{PLUS} glass columns with 10 mm and 15 mm ID.

Connection Principle:



Piston
Thread: 1/4"-28G

Ferrule

Nut


Tubing 1/16"
(0.8 x 1.6)

Connection to
Column
1/4"-28G



K.P316

Connection to
Column/System
1/4"-28G



JR-CFL-CB1PK-S

JR-55050

Connection to
System
10-32 HPLC



JR-58000

Connection to
System
M6




K.P200N

K.P207

Accessories-1/8

Accessories for ECO and ECO^{PLUS} glass columns with 25 mm – 80 mm ID.

Connection Principle:



Piston
Thread: 1/4"-28G

Ferrule

Nut


Tubing 1/8"
(1.6 x 3.2)

Connection to
Column
1/4"-28G



K.P316

Connection to
Column/System
1/4"-28G



JR-CFL-CB2PK-S

JR-55051

Connection to
System
M6



K.P300N

K.P307