




Easy Connection of YMC Glass Columns to Äkta systems

To allow easy connection of YMC glass columns to Äkta systems, YMC have specifically designed an adapter. This adapter will allow direct connection of the 10 – 32 fittings of the Äkta systems to the inlets and outlets of YMC glass columns.

This adapter is compatible to all existing YMC glass columns with internal diameters between 10 to 50 mm and can be ordered as a pack of 2 pieces.



Product Description:

- Adapter 10 – 32 on ¼ inch - 28G (connection for 10 – 32 fitting from Äkta systems), pack of 2 pieces
- Product Code: JR-0636-PCTFE/P2

For more information, please contact your YMC representative or YMC Europe GmbH

