

NEW

YMC-Actus SP

Spring-Assisted Axial Compression Columns

YMC
EUROPE GMBH

Prepacked Preparative HPLC Columns with Axial Compression Technology

- Ready-to-Use spring-assisted HPLC columns
- For all YMC packing materials with 7 µm particle size or larger
- Ideal solution for hazardous compounds
- Space-saving design

Specifications

Available IDs	50 mm (tripod stand) 80 mm (tripod or wheeled stand) 100 mm (wheeled stand)
Bed Length	up to 250 mm
Backpressure Limit	10 MPa
Compression Technology	spring-assisted axial compression
Wetted Materials	SUS316, SUS316L, PTFE



YMC-Actus SP columns are the perfect alternative to commonly used DAC columns. Hazardous compounds such as highly potent, toxic, or radioactive molecules often require dedicated systems and hardware. Prepacked YMC-Actus SP columns are the ideal solution since cleaning and validation for re-use are eliminated.

Their space-saving and wheeled design allows for maximum flexibility even with limited lab space. Another advantage: you no longer need to spend time on column packing!

YMC Europe GmbH

Schöttmannshof 19
D-46539 Dinslaken
Germany
Phone +49 2064 427-0, Fax +49 2064 427-222
www.ymc.eu

YMC Schweiz GmbH

Im Wasenboden 8
4056 Basel
Schweiz
Phone +41 61 56180-50, Fax + 41 61 56180-59
www.ymc-schweiz.ch

YMC CO., LTD.

YMC Karasuma-Gojo Bld. 284 Daigo-cho,
Karasuma Nishiiru Gojo-dori Shimogyo-ku,
Kyoto 600-8106 Japan
Phone +81 75 34245-15, Fax +81 75 34245-50
www.ymc.co.jp